

Face Shield Assembly

Each face shield is comprised of:

1. PETG face shield with an attached elastic strap
2. Foam

Assembly Instructions:

1. Remove the protective film from the shield.
2. Remove the backing of the foam to reveal the adhesive. Stick the adhesive side of the foam onto the back of the shield in line with the elastic strap.

Cleaning tips:

- Clean with a warm, gentle soap and rinse with clean water.
- A microfiber cloth is best.
- Use hydrogen peroxide (3%-5%) or isopropyl alcohol (70%).

